

 	Allegato STOOL SPECIMENT COLLECTION FOR CHEMICAL-PHYSICAL EXAMINATION	All-SDF.LAB03.38 <hr/> Pagina 1 di 1 <hr/> Aggiornato il: 05/07/2024
--	--	---

It is recommended to avoid using glycerin suppositories within 24 hours prior to collection.

The patient should defecate in a clean container.

At least 3-5 g of faeces (size of a walnut) or 5-10 ml (in case of liquid faeces) are required. Stool specimens should be placed with the appropriate small paddle into the containers.

The stool sample should not be contaminated with urine.

Physical chemical examination: stool should be delivered to the laboratory within 2 hours from the collection or stored for 24 hours in a refrigerator.