

  LABORATORIO ANALISI	Allegato URINE: ALKALIZED 24-HOUR URINE COLLECTION	All-SDF.LAB03.44 Pagina 1 di 1 Aggiornato il: 05/07/2024
---	---	---

Collect the container with the preservative from our laboratory.

- Alkalinization procedure is carried out before starting the collection by pouring 4 tablespoons of baking soda into a single graduated, well-cleaned container.

- Discard the urine from the first morning urination and write down the time.

- Collect the urine from all subsequent urinations in a single, well-cleaned graduated container
Keep the container in a cool place throughout the collection.

- The next morning complete the collection with the urine produced at the same time of the previous day.

- When ending the collection, measure the total quantity of the collected urine and thoroughly mix, pour the urine sample into the appropriate standard container, which can be purchased at the pharmacy, writing down the total volume measured.