

  LABORATORIO ANALISI	Allegato Stool Speciment Collection For Coproculture and Other Tests	All-SDF.LAB03.52 Pagina 1 di 1 Aggiornato il: 05/07/2024
---	---	---

**COPROCULTURE – HELYCOBACTER PYLORI, FAECAL ANTIGEN -
CLOSTRIDIUM DIFFICILE –MYCETES – FAECAL CALPROTECTIN**

INSTRUCTIONS FOR THE PATIENT

- Provide a container for stool collection.
- Write down the patient’s last name, first name, date of birth.
- If possible, carry out the collection in the morning of the day the specimen will be delivered to the laboratory or the previous afternoon, taking care to store the sample in a refrigerator for 24 hours.
- Dispose of the stool on a dry and clean surface.
- Collect with the disposable spoon or stick of the container a small amount of stool as large as a nut in case of formed stool, or one third of the volume of the container itself in case of liquid or broken stool (the container must not be completely filled).
- Close the container carefully, deliver to the Sampling Centre/Point or directly to the Analysis Laboratory.